



AUGENCHIRURGIE


*Für Qualität
&
Sicherheit*






www.sutures.be


 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/2 circle 	2 x HSP 8 mm 200 µm		0.5	7/0	45 cm	<i>violet</i>	OSU07
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i>	OSU061


 Side cutting spatula double armed	Metric	USP	Length	Order number			
3/8 circle 	2 x DSP 5.5 mm 150 µm		0.2	10/0	30 cm	<i>violet</i>	OSU102
	2 x DSP 6 mm 200 µm		0.4	8/0	30 cm	<i>violet</i>	OSU08-30
			0.4	8/0	45 cm	<i>violet</i>	OSU08
	2 x DSP 6.2 mm 150 µm		0.2	10/0	30 cm	<i>violet</i>	OSU10
	2 x DSP 6.4 mm 200 µm		0.3	9/0	30 cm	<i>violet</i>	OSU09
			0.5	7/0	45 cm	<i>violet</i>	OSU072
	2 x DSP 8 mm 200 µm		0.5	7/0	45 cm	<i>violet</i>	OSU071
2 x DSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i>	OSU062	
2 x DSP 10 mm 200 µm		0.5	7/0	45 cm	<i>violet</i>	OSU073	


 Side cutting spatula	Metric	USP	Length	Order number			
3/8 circle 	DSP 6.2 mm 150 µm		0.2	10/0	10 cm	<i>violet</i>	OSU101
	DSP 6.4 mm 200 µm		0.3	9/0	30 cm	<i>violet</i>	OSU09-1


 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/4 circle 	2 x VSP 8 mm 250 µm		0.7	6/0	45 cm	<i>violet</i>	OSU06-250
	2 x VSP 8 mm 300 µm 350 µm		0.7	6/0	45 cm	<i>violet</i>	OSU06
			1	5/0	45 cm	<i>violet</i>	OSU05


 Micro needle taper point double armed			Metric	USP	Length		Order number
3/8 circle 	MDR 6.4 mm		0.3	9/0	30 cm	<i>violet</i>	OSU095
	150µm						
 Reverse cutting			Metric	USP	Length		Order number
3/8 circle 	DS 12 mm		0.7	6/0	45 cm	<i>violet</i>	11071512
			1	5/0	45 cm	<i>violet</i>	11101512
	DS 12 mm		0.7	6/0	50 cm	<i>undyed</i>	11071512UP
	Easy Pass		1	5/0	50 cm	<i>undyed</i>	11101512UP


 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/2 circle 	2 x HSP 8 mm 200 µm	 	0.4	8/0	45 cm	<i>violet</i>	O910081
	2 x HSP 8 mm 200 µm	 	0.5	7/0	45 cm	<i>violet</i>	O91007
	2 x HSP 8 mm 300 µm	 	0.7	6/0	45 cm	<i>violet</i>	O910061
	2 x HSP 8 mm 300 µm	 	0.7	6/0	45 cm	<i>undyed</i>	O910061U

 Side cutting spatula double armed	Metric	USP	Length	Order number			
3/8 circle 	2 x DSP 6 mm 200 µm	 	0.4	8/0	45 cm	<i>violet</i>	O91008
	2 x DSP 6 mm 220 µm	 	0.5	7/0	45 cm	<i>violet</i>	O910072
	2 x DSP 6.5 mm 200 µm	 	0.4	8/0	30 cm	<i>violet</i>	O910082
	2 x DSP 6.5 mm 200 µm	 	0.5	7/0	30 cm	<i>violet</i>	O910077
	2 x DSP 11 mm 350 µm	 	1	5/0	45 cm	<i>violet</i>	O910051


 Side cutting spatula	Metric	USP	Length	Order number			
3/8 circle 	DSP 6.5 mm 200 µm		0.4	8/0	30 cm	<i>violet</i>	O910082-1
	DSP 6.5 mm 200 µm		0.5	7/0	30 cm	<i>violet</i>	O910077-1


 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/4 circle 	2 x VSP 6.3 mm 250 µm	 	0.7	6/0	45 cm	<i>violet</i>	O910063
	2 x VSP 8 mm 300 µm	 	0.7	6/0	45 cm	<i>violet</i>	O91006
	2 x VSP 8 mm 300 µm	 	0.7	6/0	20 cm	<i>undyed</i>	O91006U-20
	2 x VSP 8 mm 300 µm	 	0.7	6/0	45 cm	<i>undyed</i>	O91006U
	2 x VSP 8 mm 350 µm	 	1	5/0	45 cm	<i>violet</i>	O91005
	2 x VSP 8 mm 350 µm	 	1	5/0	45 cm	<i>undyed</i>	O91005U

 Side cutting spatula double armed	Metric	USP	Length	Order number			
3/8 circle 	2 x DSP 6.5 mm 200 µm		0.4	8/0	45 cm	<i>black</i>	OS81B
	2 x DSP 6.5 mm 240 µm		0.5	7/0	45 cm	<i>black</i>	OS71
			0.7	6/0	45 cm	<i>black</i>	OS61
2 x DSP 8 mm 350 µm		1	5/0	30 cm	<i>black</i>	OS52	


 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/4 circle 	2 x VSP 8 mm 250 µm		0.7	6/0	45 cm	<i>black</i>	OS62
	2 x VSP 10 mm 350 µm		1	5/0	45 cm	<i>black</i>	OS51

 Reverse cutting	Metric	USP	Length	Order number			
3/8 circle 	DS 12 mm		0.7	6/0	45 cm	<i>black</i>	8071512
			1	5/0	45 cm	<i>black</i>	8101512
			1.5	4/0	75 cm	<i>black</i>	8151512
	DS 16 mm		1.5	4/0	75 cm	<i>black</i>	8151516
	DS 19 mm		1.5	4/0	75 cm	<i>black</i>	8151519


 Reverse cutting double armed	Metric	USP	Length	Order number			
3/8 circle 	2 x DS 12 mm		0.7	6/0	45 cm	<i>black</i>	8071512-2

 Side cutting spatula double armed	Metric	USP	Length	Order number
3/8 circle  2 x DSP 6.2 mm 150 µm 	0.4	8/0	45 cm	<i>blue</i> OS82
2 x DSP 6.5 mm 200 µm 	0.4	8/0	45 cm	<i>blue</i> OS81

 **Round bodied taper point double armed** **Metric** **USP** **Length** **Order number**

1/4 circle 	2 x VR 8 mm 220 µm		0.4	8/0	45 cm	<i>black</i>	ON81
--	-----------------------	---	-----	-----	-------	--------------	------


 **Side cutting spatula** **Metric** **USP** **Length** **Order number**


1/2 circle 	HSP 5.5 mm 150 µm		0.2	10/0	30 cm	<i>black</i>	ON1081
--	----------------------	---	-----	------	-------	--------------	--------


 **Side cutting spatula double armed** **Metric** **USP** **Length** **Order number**


1/2 circle 	2 x HSP 5.5 mm 150 µm		0.2	10/0	30 cm	<i>black</i>	ON108
	2 x HSP 6.2 mm 150 µm		0.2	10/0	30 cm	<i>black</i>	ON102
	2 x HSP 6.4 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON106
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>black</i>	ON62



 **Side cutting spatula** **Metric** **USP** **Length** **Order number**


3/8 circle 	DSP 4.75 mm 150 µm		0.2	10/0	13 cm	<i>black</i>	ON1021
			0.3	9/0	13 cm	<i>black</i>	ON99
	DSP 6.2 mm 150 µm		0.2	10/0	15 cm	<i>black</i>	ON1011


 Side cutting spatula double armed	Metric	USP	Length	Order number			
3/8 circle 	2 x DSP 6 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON103
			0.1	11/0	20 cm	<i>black</i>	ON111
			0.2	10/0	30 cm	<i>black</i>	ON101
	2 x DSP 6.2 mm 150 µm		0.3	9/0	30 cm	<i>black</i>	ON95
			0.3	9/0	45 cm	<i>black</i>	ON95-45
			0.4	8/0	30 cm	<i>black</i>	ON85
			0.4	8/0	45 cm	<i>black</i>	ON85-45
	2 x DSP 6.2 mm 200 µm		0.5	7/0	30 cm	<i>black</i>	ON75
	2 x DSP 6.4 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON105
			0.2	10/0	38 cm	<i>black</i>	ON105-38
		0.2	10/0	45 cm	<i>black</i>	ON105-45	
2 x DSP 8 mm 300 µm		0.7	6/0	45 cm	<i>black</i>	ON61	


 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/4 circle 	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>black</i>	ON51


 Reverse cutting	Metric	USP	Length	Order number			
3/8 circle 	DS 12 mm		0.5	7/0	75 cm	<i>black</i>	9051512
			0.7	6/0	75 cm	<i>black</i>	9071512


 Side cutting spatula double armed	Metric	USP	Length	Order number	
1/2 circle  2 x HSP 8 mm 350 µm 	1	5/0	45 cm	<i>white</i> OP52	
3/8 circle  2 x DSP 8 mm 350 µm 	1	5/0	45 cm	<i>white</i> OP53	
1/4 circle 	2 x VSP 8 mm 250 µm 	0.5	7/0	45 cm	<i>white</i> OP71
	2 x VSP 8 mm 300 µm 	0.7	6/0	45 cm	<i>white</i> OP61
	2 x VSP 8 mm 350 µm 	1	5/0	45 cm	<i>white</i> OP51
	2 x VSP 8 mm 350 µm 	1.5	4/0	45 cm	<i>white</i> OP41



 Side cutting spatula double armed	Metric	USP	Length	Order number
3/8 circle  2 x DSP 6.4 mm 200 µm	0.2	10/0	30 cm	<i>blue</i> OPP105
	0.3	9/0	30 cm	<i>blue</i> OPP95

 Side cutting spatula	Metric	USP	Length	Order number
3/8 circle  DSP 6 mm 200 µm	0.5	7/0	45 cm	<i>blue</i> OPP71-1

 Side cutting spatula double armed	Metric	USP	Length	Order number
1/4 circle  2 x VSP 8 mm 300 µm	0.7	6/0	45 cm	<i>blue</i> OPP62
	0.2	10/0	20 cm	<i>blue</i> OPP116

 Side cutting spatula	Metric	USP	Length	Order number
1/4 circle  VSP 16 mm 150 µm	0.2	10/0	10 cm	<i>blue</i> OPP103

 Side cutting spatula loop double armed	Metric	USP	Length	Order number
1/4 circle  VSP 15.3 mm 150 µm	0.2	10/0	20 cm	<i>blue</i> OPP113L-2
	0.2	10/0	20 cm	<i>blue</i> OPP103L-2

 Side cutting spatula loop	Metric	USP	Length	Order number
1/4 circle  VSP 16 mm 150 µm 	0.2	10/0	20 cm	<i>blue</i> OPP103L
 Side cutting spatula double armed	Metric	USP	Length	Order number
Straight  2 x GSP 16.15 mm 150 µm 	0.2	10/0	20 cm	<i>blue</i> OPP101
	0.3	9/0	20 cm	<i>blue</i> OPP91
 Side cutting spatula loop	Metric	USP	Length	Order number
Straight  GSP 16.15 mm 150 µm 	0.2	10/0	15 cm	<i>blue</i> OPP110
 Side cutting spatula double armed	Metric	USP	Length	Order number
Straight  GSP 16 mm 150 µm DSP 4 mm 100 µm 	0.2	10/0	10 cm	<i>blue</i> OPP102
3/8 circle  GSP 16 mm 150 µm DSP 5 mm 150 µm 	0.2	10/0	20 cm	<i>blue</i> OPP115